

Anabolic Lab™

Analytical Lab Testing for
Anabolic Steroid Harm Reduction

<https://anaboliclab.com>

Lab Report – April 27, 2016

Balkan Pharmaceuticals
Nandrolone D (200mg/mL)
Nandrolone Decanoate



If you like the attached report, please visit [our website](https://anaboliclab.com) to find out how you can support independent and objective lab testing.



Simec AG

Areal Bleiche West Postfach 413 CH - 4800 Zofingen
Telefon +41 62 752 83 08 Fax +41 62 752 83 09
E-Mail: info@simec.ch www.simec.ch

Test report

Number **C 42756**

Object of analysis: Balkan Pharma Nandrolone D

Customer: AnabolicLab.com

Samples: 1

Batch: ND003C, Exp. 02/2019

Receipt: 01.04.2016


Subject: Assay, Identification Nandrolone Decanoate

Method(s): HPLC

Sample / Analysis:	Method	Status	Result	Specification	complies
<u>Nandrolone Decanoate:</u>					
Identification	HPLC	N	complies	complies	yes
Content	HPLC	N	197.6 mg/mL	200 mg/mL	-

Method-status: G = GMP A = accredited V = generally validated P = validated on product N = not validated
E = external lab

Remarks: The sample(s) mentioned above have been analysed as they have been sent to us by the client. We have not controlled external processes like manufacturing, labelling, sampling, shipping, storage. These results are for information only and do not compensate for correct quality control by the manufacturer/distributor. Generally, pharmaceutical products have to be produced and distributed under full GMP/GDP regime, including GMP-compliant analyses with validated methods.

Signatures:  **controlled & released:**

Staff: 

completed:


27.04.2016

The enclosed results refer exclusively to the object of examination described above. The measuring accuracy is available on request. Only handwritten signatures are valid. Without written permission, it is not allowed to publish single parts of this report.

page 1 of 1

Report steroid analysis

HPLC: Nr.5
 Analysendatum: 22.04.2016 15:09:34 Druckdatum: 25.04.2016 08:56:02
 Vial: 25 Injektvolumen: 20
 Probenidentifikation: C42756 A
 Operator: admin
 Chrom.-Beding. ACN / Wasser, gradient
 Temperatur: 40 °C
 Säule Nr. 121, 250 x 4,6, C18
 Detektion: UV 240nm



Channel A

UV Results

Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	40,840	93839784	0,000	100,000


Totals		93839784	0,000	100,000
--------	--	----------	-------	---------

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2016\April 2016\Diverse Steroide\19,04,16\032_C42756 A-Rep1.dat
 Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
 Analysendatum: 22.04.2016 16:02:53
 Vial: 25
 Probenidentifikation: C42756 A
 Operator: admin
 Chrom.-Beding. ACN / Wasser, gradient
 Temperatur: 40 °C
 Säule Nr. 121, 250 x 4,6, C18
 Detektion: UV 240nm

Druckdatum: 25.04.2016 08:56:06
 Injektvolumen: 20




Channel A

UV Results				
Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	40,897	93882940	0,000	98,405
Totals		95405067	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2016\April 2016\Diverse Steroide\19,04,16\032_C42756 A-Rep2.dat
 Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
 Analysendatum: 22.04.2016 16:56:16 Druckdatum: 25.04.2016 08:56:10
 Vial: 25 Injektvolumen: 20
 Probenidentifikation: C42756 A
 Operator: admin
 Chrom.-Beding. ACN / Wasser, gradient
 Temperatur: 40 °C
 Säule Nr. 121, 250 x 4,6, C18
 Detektion: UV 240nm




Channel A

UV Results				
Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	41,073	93303138	0,000	98,700
Totals		94531630	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2016\April 2016\Diverse Steroide\19,04,16\032_C42756 A-Rep3.dat
 Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC:	Nr.5	Druckdatum:	25.04.2016 08:56:14
Analysendatum:	22.04.2016 17:49:40	Injektvolumen:	20
Vial:	26		
Probenidentifikation:	C42756 B		
Operator:	admin		
Chrom.-Beding.	ACN / Wasser, gradient		
Temperatur:	40 °C		
Säule	Nr. 121, 250 x 4,6, C18		
Detektion:	UV 240nm		



Channel A

UV Results


Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	41,063	91679201	0,000	98,373
Totals		93195487	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2016\April 2016\Diverse Steroide\19,04,16\033_C42756 B-Rep1.dat
 Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
 Analysendatum: 22.04.2016 18:43:07
 Vial: 26
 Probenidentifikation: C42756 B
 Operator: admin
 Chrom.-Beding. ACN / Wasser, gradient
 Temperatur: 40 °C
 Säule Nr. 121, 250 x 4,6, C18
 Detektion: UV 240nm

Druckdatum: 25.04.2016 08:56:18
 Injektvolumen: 20



Channel A

UV Results

Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	41,073	91544909	0,000	98,400


Totals		93033567	0,000	100,000
--------	--	----------	-------	---------

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2016\April 2016\Diverse Steroide\19,04,16\033_C42756 B-Rep2.dat
 Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC: Nr.5
 Analysendatum: 22.04.2016 19:36:31
 Vial: 26
 Probenidentifikation: C42756 B
 Operator: admin
 Chrom.-Beding. ACN / Wasser, gradient
 Temperatur: 40 °C
 Säule Nr. 121, 250 x 4,6, C18
 Detektion: UV 240nm

Druckdatum: 25.04.2016 08:56:22
 Injektvolumen: 20




Channel A

UV Results				
Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	41,103	91603774	0,000	98,344
Totals		93146090	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2016\April 2016\Diverse Steroide\19,04,16\033_C42756 B-Rep3.dat
 Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met

Report steroid analysis

HPLC:	Nr.5	Druckdatum:	25.04.2016 08:56:26
Analysendatum:	22.04.2016 20:29:52	Injektvolumen:	20
Vial:	1		
Probenidentifikation:	Blank		
Operator:	admin		
Chrom.-Beding.	ACN / Wasser, gradient		
Temperatur:	40 °C		
Säule	Nr. 121, 250 x 4,6, C18		
Detektion:	UV 240nm		



Channel A

UV Results				
Name	Time	Area	Amount	Area Percent
Nandrolone decanoate	16,247	1514010	0,000 0,000 BDL	100,000
Totals		1514010	0,000	100,000

Daten File: C:\EZChrom Elite\Enterprise\Projects\Simec\Data\2016\April 2016\Diverse Steroide\19,04,16\034_Blank.dat
 Methode: C:\EZChrom Elite\Enterprise\Projects\Simec\Method\Steroide\Steroide A.met